



The CryoSafe CryoMizer Series Tanks

Models	CM-1	CM-2	CM-3	CM-4
Vial Capacity	875	2,500	5,000	7,000
Box Capacity	35	25	50	70
Rack Capacity	7	5	5	7
Vials per Box	25	100	100	100
LN2 Capacity (Liters)	47.4	61	121	175
Static Evaporation Rate (Liters/Day)*	0.39	0.85	0.99	0.99
Normal Working Time (Days)**	76	38	70	104
Neck Diameter in (mm)	5 (127)	8.5 (216)	8.5 (216)	8.5 (216)
Exterior Height in (mm)	26.5 (673)	27.25 (692)	38 (965)	37.75 (959)
Exterior Diameter in (mm)	20 (508)	22 (559)	22 (559)	26 (665)