



## The CryoSafe Liquid Dewar Series

Models	LD 4	LD 5	LD 10	LD 20	LD 30	LD 50
LN2 Capacity (Liters)	4	5	10	20	32	50
Static Evaporation Rate (Liters/Day)	0.19	0.15	0.18	0.18	0.22	0.49
Neck Diameter in (mm)	1.4 (35.5)	2.18 (56)	2.18 (56)	2.18 (56)	2.5 (64)	2.5 (64)
Exterior Height in (mm)	16.8 (426)	18.2 (462)	21.5 (546)	24.5 (622)	24 (610)	30.5 (775)
Exterior Diameter in (mm)	7.3 (185)	8.8 (222)	10.3 (260)	14.5 (368)	17 (432)	17 (432)
Weight Empty lbs (kg)	6 (2.7)	8 (4)	12 (5.4)	19 (9)	25 (11.4)	31 (14)
Weight Full lbs (kg)	13 (6)	17 (8)	31 (14)	55 (25)	82 (37.2)	120 (54.4)