

# L A Y O U T S

## Auto-Fill Inventory Systems

Rack / Box

Frame / Cassette



### 21" DIAMETER TANKS

7 each 13-2 Tower Racks  
91 each 100/81 Cell Boxes  
4 each 13-2 Mini Tower Racks  
52 each 25 Cell Boxes  
Total Vial Capacity 10,400 (2ml)

56 each Frames\*  
224 each Cassette\*



### 25" DIAMETER TANKS

12 each 13-2 Tower Racks  
156 each 100/81 Cell Boxes  
4 each 13-2 Mini Tower Racks  
52 each 25 Cell Boxes  
Total Vial Capacity 16,900 (2ml)

82 each Frames\*  
328 each Cassette\*



### 31" DIAMETER TANKS

17 each 13-2 Tower Racks  
221 each 100/81 Cell Boxes  
6 each 13-2 Mini Tower Racks  
78 each 25 Cell Boxes  
Total Vial Capacity 24,050 (2ml)


140 each Frames\*  
560 each Cassette\*



### 32" DIAMETER TANKS

18 each 13-2 Tower Racks  
234 each 100/81 Cell Boxes  
6 each 13-2 Mini Tower Racks  
78 each 25 Cell Boxes  
Total Vial Capacity 25,350 (2ml)


144 each Frames\*  
576 each Cassette\*



### 39" DIAMETER TANKS

28 each 13-2 Tower Racks  
364 each 100/81 Cell Boxes  
6 each 13-2 Mini Tower Racks  
78 each 25 Cell Boxes  
Total Vial Capacity 38,350 (2ml)


228 each Frames\*  
912 each Cassette\*



### CryoGuard 1

14 each 13-2 Tower Racks  
182 each 100/81 Cell Boxes  
4 each 13-2 Mini Tower Racks  
52 each 25 Cell Boxes  
Total Vial Capacity 19,500 (2ml)


80 each Frames\*  
400 each Cassette\*



### CryoGuard 2

30 each 13-2 Tower Racks  
390 each 100/81 Cell Boxes  
8 each 13-2 Mini Tower Racks  
104 each 25 Cell Boxes  
Total Vial Capacity 41,600 (2ml)

164 each Frames\*  
820 each Cassette\*



### CryoGuard 3

58 each 13-2 Tower Racks  
754 each 100/81 Cell Boxes  
12 each 13-2 Mini Tower Racks  
156 each 25 Cell Boxes  
Total Vial Capacity 79,300 (2ml)

301 each Frames\*  
1505 each Cassette\*



\*Fenwal Blood Bag 4R9953