

Product Name:	CM 3 CS
Description:	<p>CryoMizer Class Cryogenic Freezer  CryoMizer Cardboard Package System Includes:  (5) Stainless Steel Racks, Including Center Max Rack  (50) 81 or 100 Cell Cardboard Boxes and Grids  (Please specify when ordering)  Rollerbase for CM-3  One Pair CryoGloves  Liquid Level Measuring Ruler  Model Specifications:  5,000 (2ml) Total Vial Capacity  50 Box Capacity (100 or 81 Vials per box) (Boxes not included)  5 Rack Capacity (10 boxes per rack)  121 Liter Liquid Nitrogen Capacity  Static Evaporation Rate (Liters/Days): 0.99  Normal Working Time: 70 days  Static Holding Time: 122 days  Neck Diameter (in/mm): 8.5/215.9  Exterior Diameter (in/mm): 22/558.8  Exterior Height (in/mm): 38/965.2  Weight Empty (lbs/kg): 81/36.7  Weight Full (lbs/kg): 300/136.1  Shipping Weight (lbs/kg) includes rollerbase: 200/90.7  Three Year Parts and Labor Warranty  Five Years Vacuum Warranty</p>