

Product Name:	LD 10
Description:	CryoSafe Liquid Nitrogen Dewar Model specifications: Liquid nitrogen capacity (Liters): 10 Static evaporation rate (Liters/Day): 0.18 Neck diameter (in/mm): 2.18/55.4 Exterior height (in/mm): 21.5/546.1 Exterior diameter (in/mm): 10.3/261.6 Weight full (lbs/kg): 31/14.1 Weight empty (lbs/kg): 12/5.4 Three year parts and labor warranty Five year vacuum warranty